



Apnea Monitor

Model #: _____
Serial #: _____
Date: _____

Provider Name: _____
Phone #: _____
Contact Person: _____

Setting up your equipment:

- Place the monitor where it can be heard throughout your home.
- Plug the unit into a 3-prong grounded outlet only.
- Make sure to have plenty of lead wire to avoid pulling equipment.

Using your equipment:

- The settings on the monitor will be pre-set by a respiratory therapist as prescribed by your physician.
- Place electrodes under each arm or as shown by therapist.
- Connect lead wires to the electrode then plug into the monitor. The alarms will be explained at the time of delivery.

Maintenance:

- Keep your monitor free from oxygen sources, liquids, dirt, and other debris.
- Have your provider test the monitor every month to verify it is in proper working order and you have adequate supplies.
- Skin needs to be clean and dry to allow for proper electrode contact. This will also decrease frequency of false alarms. Keep skin clean by using a mild cleansing agent.

Safety Issues:

- Be sure you know proper procedure in case of alarm signal.
- Respond to all alarms.
- Place your monitor where it will be heard.
- Check for frayed wires and possible shock sources.
- Have emergency procedure posted clearly along with the phone number of your provider.
- Be sure the battery charger is plugged into the wall outlet securely to ensure a good charge.
- Do not place monitor in or near water.
- Check lead wires to avoid entanglement (during infant use).

Frequent replacement items:

- Electrodes- monthly to bimonthly
- Wires- monthly
- Belts- as needed

Please contact your equipment provider for any questions, adjustments, or repairs.

Your local provider:



©Copyright The VGM Group 2002
B.L.U.E. Sheet photos, drawings, and text courtesy of
Golden Technologies, Invacare, Leisure-Lift, Permobil,
Puritan Bennett, Pride, Respironics, Ross,
Sammons Preston, Smith-Nephew, Sunrise Medical