



Commodes

Model #: _____

Provider Name: _____

Serial #: _____

Phone #: _____

Date: _____

Contact Person: _____

Setting up your equipment:

- Make sure commode is adjusted to proper height for user.

Using your equipment:

- Seat must be down before using commode; do not sit on pail holder. (Diagram 1)
- If using a mobile commode, be sure wheels are locked before using. (Diagram 2)
- If commode has drop-arms, ensure they are locked in the upright position before using. (Diagram 3)
- If using elevated commode seat on toilet bowl, make sure the seat is properly attached to prevent loss of balance and/or tipping.

Diagram 1

Pail
Holder



Diagram 2

Wheels
with locks



Diagram 3

Drop Arm
Commode



Frequent replacement items:

- Tips on leg extensions

Maintenance:

- Make sure all screws; nuts and bolts are tight at all times.
- Inspect rubber tips on the legs to ensure they are not worn or torn.
- For hygiene sake, empty pail and clean it after each use.
- Make sure wheels and locks are working properly.

Safety Issues:

- Do not sit on commode for extended periods of time.
- Commode seats are not intended to serve as wheelchairs.
- Make sure the user can support themselves while using the commode. If user is unsteady or weak, a caregiver should provide support during the toileting process.

Please contact your equipment provider for any questions, adjustments, or repairs.

Your local provider:



©Copyright The VGM Group 2002
B.L.U.E. Sheet photos, drawings, and text courtesy of
Golden Technologies, Invacare, Leisure-Lift, Permobil,
Puritan Bennett, Pride, Respironics, Ross,
Sammons Preston, Smith-Nephew, Sunrise Medical