

Model #: _____
Serial #: _____
Date: _____

Provider Name: _____
Phone #: _____
Contact Person: _____

Things to Remember While Using Oxygen:

- Change nasal cannula and/or oxygen mask every two weeks. (Picture 1)
- Change oxygen tubing every three months. (Picture 2)
- Clean humidifier bottle twice weekly with hot soapy water. If you have a disposable humidifier bottle, it should be changed monthly. (Picture 3)
- Call your physician if you experience frequent headaches, anxiety, blue lips or fingernails, drowsiness, confusion, restlessness, or slow, shallow, or irregular breathing.

Picture 1



Picture 2



Picture 3



Safety Issues:

- DO NOT smoke while using oxygen.
- DO NOT stand within five feet of someone who is smoking while you are on oxygen.
- Keep NO SMOKING signs posted in your place of residence.
- Keep heat sources at least five feet away from you and your oxygen tubing.
- DO NOT use an oil-based product (i.e. Petroleum jelly) for nose irritation, instead use a water-based product.
- Keep oxygen cylinders stored in a well-ventilated area. If cylinders cannot be secured in an upright position, lay them on their sides.
- DO NOT change your oxygen flow rate without consulting your physician.

Please contact your equipment provider for any questions, adjustments, or repairs.

Your local provider:



©Copyright The VGM Group 2002

B.L.U.E. Sheet photos, drawings, and text courtesy of
Golden Technologies, Invacare, Leisure-Lift, Permobil,
Puritan Bennett, Pride, Respironics, Ross,
Sammons Preston, Smith-Nephew, Sunrise Medical