



TENS Therapy

Model #: _____
Serial #: _____
Date: _____

Provider Name: _____
Phone #: _____
Contact Person: _____

Setting up your equipment:

- Install the battery.
- Connect one pair of electrodes to each output channel using the lead wires supplied.



Using your equipment:

1. Clean the skin where the electrodes are to be placed.
2. If the electrodes do not contain gel, then gel should be applied to the skin where the electrodes will be applied.
3. Place the electrodes on the desired area according to your physician's instructions.
4. Turn the unit on.
5. Adjust the intensity, duration, and number of pulses as directed by your physician.
6. Turn unit off before removing the electrodes from your skin.

Frequent Replacement Items:

- Lead wires (every 6 months or as needed)
- Electrodes (every month)
- Disposable batteries



Maintenance:

- Clean the case with isopropyl alcohol or cleaning agent.
- Do not submerge the device in liquid or use excessive liquid to clean the unit.

Safety Issues:

- TENS devices can adversely affect the operation of demand-type cardiac pacemakers.
- TENS is not recommended for patients with known heart disease without a physician's evaluation of risk.
- Do not stimulate over the eyes or carotid sinus nerves.
- Do not place electrodes in a manner that causes current to flow through the head.
- Electronic equipment such as ECG monitors and ECG alarms may not operate properly when TENS is in use.
- Avoid adjusting controls while operating machinery or vehicles.



Please contact your equipment provider for any questions, adjustments, or repairs.

Your Local Provider:



©Copyright The VGM Group 2005

B.L.U.E. Sheet photos, drawings, and text courtesy of BioMedical Life, Golden Technologies, Invacare, Leisure-Lift, Permobil, Puritan Bennett, Pride, Rehabicare, Respironics, Ross, Sammons Preston, Smith-Nephew, Sunrise Medical